









Instruction for Use

**Ambu<sup>®</sup> aScope<sup>™</sup> 4  
RhinoLaryngo Slim**

## Contents

## Page

English (Instruction for Use).....	4
Български (Инструкция за употреба).....	11
Čeština (návod k použití).....	18
Dansk (Brugsvejledning).....	25
Deutsch (Bedienungsanleitung).....	32
Ελληνικά (Οδηγίες χρήσης).....	39
Español (Instrucciones de uso).....	46
Eesti (Kasutusjuhend).....	53
Suomi (Käyttöohjeet).....	60
Français (Mode d'emploi).....	67
Hrvatski (Upute za upotrebu).....	74
Magyar (Használati útmutató).....	81
Italiano (Istruzioni per l'uso).....	88
英語 (使用説明書).....	95
Lietuvių (naudojimo instrukcija).....	102
Latviešu val. (Lietošanas instrukcija).....	109
Nederlands (Gebruiksaanwijzing).....	116
Norsk (Bruksanvisning).....	123
Polski (Instrukcja użytkowania).....	130
Português (Instruções de utilização).....	137
Română (Instrucțiuni de utilizare).....	144
Русский (Инструкция по применению).....	151
Slovensky (Návod na použitie).....	158
Slovenščina (Navodila za uporabo).....	165
Svenska (bruksanvisning).....	172
Türkçe (Kullanma Talimatı).....	179
中文 (使用说明).....	186

**CE** 0086  
Pat. Pending

## 1. Important Information – Read Before Use

Read the safety instructions carefully before using the Ambu® aScope™ 4 RhinoLaryngo Slim. The instructions for use may be updated without further notice. Copies of the current version are available upon request. Please be aware that these instructions do not explain or discuss clinical procedures. They describe only the basic operation and precautions related to the operations of aScope 4 RhinoLaryngo Slim.

Before initial use of the aScope 4 RhinoLaryngo Slim it is essential for operators to have received sufficient training in clinical endoscopic techniques and to be familiar with the intended use, warnings and cautions mentioned in these instructions.

In this Instruction for Use, the term *endoscope* refers to instructions for aScope 4 RhinoLaryngo Slim, and *system* refers to aScope 4 RhinoLaryngo Slim and the compatible Ambu monitor. This Instruction for Use applies for the endoscope and information relevant for the system.

### 1.1. Intended Use

The endoscope is a sterile, single-use, flexible endoscope intended for endoscopic procedures and examination within the nasal lumens and upper airway anatomy. The endoscope is intended to provide visualization via a monitor.

The endoscope is intended for use in a hospital environment. It is designed for use in adults.

### 1.2. Warnings and Cautions

Failure to observe these warning and cautions may result in patient injury or damage to the equipment.

Ambu is not responsible for any damage to the system or patient resulting from incorrect use.

## WARNINGS

1. Only to be used by physicians, trained in clinical endoscopic techniques and procedures.
2. The endoscope is a single-use product and must be handled in a manner consistent with accepted medical practice for such devices in order to avoid contamination of the endoscope prior to insertion.
3. Do not soak, rinse, or sterilize this device as these procedures may leave harmful residues or cause malfunction of the device. Reuse of the endoscope can cause contamination, leading to infections.
4. Do not use the endoscope if the sterilisation barrier or its packaging is damaged.
5. Do not use the endoscope if it is damaged in any way or if the preuse check fails (see section 4.1).
6. The images must not be used as an independent diagnostic of any pathology. Physicians must interpret and substantiate any finding by other means and in the light of the patient's clinical characteristics.
7. Do not use active endoscopic accessories such as laser probes and electrosurgical equipment in conjunction with the endoscope, as this may result in patient injury or damage to endoscope.
8. The endoscope is not to be used when delivering highly flammable anaesthetic gases to the patient. This could potentially cause patient injury.
9. Patients should be adequately monitored at all times. Failure to do so may harm the patient.
10. Always make sure that the bending section is in a straight position when inserting and withdrawing the endoscope. Do not operate the bending lever and never use excessive force, as this may result in injury to the patient and/or damage to the endoscope.
11. Do not use excessive force when advancing, operating or withdrawing the endoscope as this may result in patient injury or damage to the endoscope.
12. The distal end of the endoscope may get warm due to heating from the light emission part. Avoid long periods of contact between the tip of the device and the mucosal membrane as sustained contact with the mucosal membrane may cause mucosal injury.

## CAUTIONS

1. Have a suitable backup system readily available in case a malfunction should occur.
2. Be careful not to damage the insertion cord or distal tip. Do not allow other objects or sharp devices such as needles to strike the endoscope.
3. US federal law restricts these devices for sale only by, or on the order of, a physician.


## 2. System Description

The endoscope can be connected to the compatible monitor. For information about the compatible monitor, please refer to its Instruction for Use.

### 2.1. System Parts

#### Endoscopes

#### Part numbers:

510001000 aScope 4 RhinoLaryngo Slim

aScope 4 RhinoLaryngo Slim is not available in all countries. Please contact your local sales office.

#### Product Name

#### Colour

#### Outer Diameter [mm]

aScope 4 RhinoLaryngo Slim

Purple

min 3.0; max 3.5

#### Compatible monitors

#### Part numbers:


Ambu® aView™  
(Reusable)

405002000 Model no. JANUS2-W08-R10  
(SW versions ≥ V2.02)

For aView model no., please check the backside label on aView.  
aView is not available in all countries. Please contact your local sales office.

## 2.2. Endoscope Parts

No.	Part	Function	Material
1	Handle	Suitable for left and right hand	MABS
2	Control lever	Moves the distal tip up or down in a single plane	POM
3	Insertion cord	Flexible airway insertion cord	PU
4	Bending section	Manoeuvrable part	PU
5	Distal end	Contains the camera and light source (two LEDs)	PC
3-4-5	Insertion portion	The ensemble of insertion cord, bending section and distal end	See above
6	Connector on the endoscope cable	Connects to blue socket on the monitor	PVC
7	Endoscope cable	Transmits the image signal to the monitor	PVC
8	Protection pipe	Protects the insertion cord during transport and storage. Remove before use.	PP
-	Packaging	Sterile barrier	Cardboard, Tyvek

Abbreviations: MABS (Methyl Acrylonitrile Butadiene Styrene), PU (Polyurethane), PP (polypropylene), PC (Polycarbonate), POM (Polyoxymethylen).

## 3. Explanation of Symbols Used

Symbols for the endoscope devices	Indication
	Working length of the endoscope insertion cord.
	Maximum insertion portion width (Maximum outer diameter).
	Field of view.
	Do not use if the product sterilisation barrier or its packaging is damaged.
	This product is not made with natural rubber latex.
	Temperature limitation: between 10 °C (50 °F) and 40 °C (104 °F) in operating environment.
	Humidity limitation: relative humidity between 30 and 85 % in operating environment.
	Atmospheric pressure limitation: between 80 and 109 kPa in operating environment.
	Manufacturer.
	Consult Instruction for Use.
	CE mark. The product complies with the EU Council directive concerning Medical Devices 93/42/EEC.
	Electrical Safety Type BF Applied Part.
	Use By, followed by YYYY-MM-DD.
	Sterile Product, Sterilisation by ETO.
	Single use product, do not reuse.
	Reference Number.
	Lot Number, Batch Code.
	UL Recognized Component Mark for Canada and the United States.
	Warning.

## 4. Use of the endoscope

Optimize patient position and consider applying relevant anesthetics to minimize patient discomfort.

Numbers in gray circles below refer to illustrations on page 2.

### 4.1. Preuse check of the endoscope

1. Check that the pouch seal is intact before opening. **1a**
2. Make sure to remove the protective elements from the insertion cord. **1b**
3. Check that there are no impurities or damage on the product such as rough surfaces, sharp edges or protrusions which may harm the patient. **1c**

Refer to the **Instruction for Use for the compatible monitor for preparation and inspection of the monitor.** **2**

### 4.2. Inspection of the Image

1. Plug in the endoscope into the corresponding connector on the compatible monitor. Please ensure the colours are identical and be careful to align the arrows. **3**
2. Verify that a live video image appears on the screen by pointing the distal end of the endoscope towards an object, e.g. the palm of your hand. **4**
3. Adjust the image preferences on the compatible monitor if necessary (please refer to the monitor Instruction for Use).
4. If the object cannot be seen clearly, clean the tip.

### 4.3. Preparation of the endoscope

Carefully slide the bending control lever forwards and backwards to bend the bending section as much as possible. Then slide the bending lever slowly to its neutral position. Confirm that the bending section functions smoothly and correctly and returns to a neutral position. **5**

### 4.4. Operating the endoscope

#### Holding the endoscope and manipulating the tip **6**

The handle of the endoscope can be held in either hand. The hand that is not holding the endoscope can be used to advance the insertion cord into the patient's nose or mouth. Use the thumb to move the control lever. The control lever is used to flex and extend the tip of the endoscope in the vertical plan. Moving the control lever downward will make the tip bend anteriorly (flexion). Moving it upward will make the tip bend posteriorly (extension). The insertion cord should be held as straight as possible at all times in order to secure an optimal tip bending angle.

#### Insertion of the endoscope **7**

To ensure the lowest possible friction during insertion of the endoscope the insertion cord may be lubricated with a medical grade lubricant. If the images of the endoscope becomes unclear, clean the tip. When inserting the endoscope orally, it is recommended to use a mouthpiece to protect the scope from being damaged.

#### Withdrawal of the endoscope **8**

When withdrawing the endoscope, make sure that the control lever is in the neutral position. Slowly withdraw the endoscope while watching the live image on the monitor.

### 4.5. After Use

#### Visual check **9**

Inspect the endoscope for any evidence of damage on the bending section, lens, or insertion cord. In case of corrective actions needed based on the inspection act according to local hospital procedures.

#### Final steps **10**

Disconnect the endoscope from the Ambu monitor and dispose the endoscope in accordance with local guidelines for collection of infected medical devices with electronic components.

## 5. Technical Product Specifications

### 5.1. Standards Applied

The endoscope function conforms with:

- IEC 60601-1: Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.
- IEC 60601-1-2: Medical electrical equipment – Part 1-2 General requirements for safety – Collateral standard: Electromagnetic compatibility - Requirements for test.
- IEC 60601-2-18: Medical electrical equipment – Part 2-18: Particular requirements for the safety of endoscopic equipment.
- ISO 8600-1: Optics and photonics - Medical endoscopes and endotherapy devices – Part 1: General requirements.
- ISO 10993-1: Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process.

### 5.2. Endoscope specifications

Insertion portion	aScope 4 RhinoLaryngo Slim
Bending section <sup>1</sup> [°]	130  ,130 
Insertion cord diameter [mm, (")]	3.0 (0.12)
Distal end diameter [mm, (")]	3.5 (0.14)
Maximum diameter of insertion portion [mm, (")]	3.5 (0.14)
Minimum tracheostomy tube size (ID) [mm]	6.0
Working length [mm, (")]	300 (11.8)
Storage and transportation	aScope 4 RhinoLaryngo Slim
Transportation temperature [°C, (°F)]	10 ~ 40 (50 ~ 104)
Recommended storage temperature <sup>3</sup> [°C, (°F)]	10 ~ 25 (50 ~ 77)
Relative humidity [%]	30 ~ 85
Atmospheric pressure [kPa]	80 ~ 109
Optical System	aScope 4 RhinoLaryngo Slim
Field of View [°]	85
Depth of Field [mm]	6 - 50
Illumination method	LED
Sterilisation	aScope 4 RhinoLaryngo Slim
Method of sterilisation	ETO
Operating environment	aScope 4 RhinoLaryngo Slim
Temperature [°C, (°F)]	10 ~ 40 (50 ~ 104)
Relative humidity [%]	30 ~ 85
Atmospheric pressure [kPa]	80 ~ 109
Altitude [m]	≤ 2000

1. Please be aware that the bending angle can be affected if the insertion cord is not kept straight.
2. Storage under higher temperatures may impact shelf life.

## 6. Trouble Shooting

If problems occur with the system, please use this trouble shooting guide to identify the cause and correct the error.

Problem	Possible cause	Recommended action
No live image on the screen but user interface is present on the display or the image is frozen	The endoscope is not connected to compatible monitor.	Connect an endoscope to the blue port on the monitor.
	The monitor and endoscope have communication problems.	Restart the monitor.
	The endoscope is damaged.	Replace the endoscope with a new one.
	A recorded image is shown on the monitor screen.	Return to live image on the monitor.
Low picture quality	Blood, saliva etc. on the lens (distal tip).	If the object cannot be seen clearly, clean the tip.

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**Ambu A/S**