

Ambu[®] Disposable Face Mask



Key Benefits

- Crystal clear dome for easy observation of the patient's condition
- Dome with thumb rest for an easy grip
- Soft, shaped, airfilled cuff allows a tight fit to the patient's face
- Supplied with hook ring; hook ring can easily be removed if not needed
- The single patient use concept offers minimized risk of cross infection
- All sizes are delivered individually packed, non-sterile, in packs of 20 units

Ambu[®] Disposable Face Mask

Ambu Disposable Face Mask is designed for use with manual and automatic resuscitators and ventilators. Complete line with 6 different sizes from neonate to large adult.

Available with or without checkvalve.

Ordering specifications

Item no.	Type
000 252 051	Ambu Disposable Face Mask Neonate
000 252 052	Ambu Disposable Face Mask Infant
000 252 053	Ambu Disposable Face Mask Toddler
000 252 054	Ambu Disposable Face Mask Small adult/Child
000 252 055	Ambu Disposable Face Mask Medium adult
000 252 056	Ambu Disposable Face Mask Large adult
000 252 951	Ambu Disposable Face Mask w. checkvalve Neonate
000 252 952	Ambu Disposable Face Mask w. checkvalve Infant
000 252 953	Ambu Disposable Face Mask w. checkvalve Toddler
000 252 954	Ambu Disposable Face Mask w. checkvalve Small adult/Child
000 252 955	Ambu Disposable Face Mask w. checkvalve Medium adult
000 252 956	Ambu Disposable Face Mask w. checkvalve Large adult



Ambu Disposable Facemasks with checkvalve

Materials

Product	Material
Mask dome	Polyvinylchloride
Cushion	Polyvinylchloride

Size

Type	Weight
Neonate	9 g
Infant	14 g
Toddler	25 g
Small adult/Child	33 g
Medium adult	47 g
Large adult	49 g



Ambu Disposable Facemasks connected to Ambu SPUR II

Ambu Australia Pty Ltd
Customer Services
 Tel.: 1300 233 118
 Fax. 1300 156 526
 aussales@ambu.com
 www.ambuaustralia.com.au

Ambu Asia (Sales office)
Customer Services
 Tel.: +604 252 9100
 Fax. +604 644 8932
 csasia@ambu.com