

Ambu® UltraSeal™ II

DATASHEET

A disposable face mask developed specifically for anaesthesia departments.

The Ambu® UltraSeal™ II is designed for use with manual and automatic ventilators and anaesthesia machines.

Ambu® UltraSeal™ II

The UltraSeal II has an easy grip, with a soft anatomically shaped air-filled cuff to facilitate a tight fit to the patient's face and enable optimal ventilation.

The mask is single patient use, so there is no need for time-consuming cleaning, disinfection or sterilisation.

Available in six different sizes from Neonate size 1 through to Large adult size 6.

All sizes come with a check valve allowing control of cuff inflation.



KEY BENEFITS

- Very soft anatomically shaped cuff enabling a tight seal with minimum applied pressure
- Flexible, self-expanding dome
- Clear dome for easy observation of the patient's condition
- All sizes are delivered individually packed in a transparent, easy to open bag
- Latex free
- DEHP and BPA free
- Check valve makes it possible to inflate and deflate the cuff

Ambu

PRODUCT SPECIFICATIONS

MATERIALS

Mask dome	Polyvinylchloride (PVC), DEHP-free
Cushion	PVC, DEHP-Free
Valve	PVC, Polyester, Nitrile, Stainless Steel

MATERIALS PACKAGING

Pouch	PE (Polyethylene)
Box	Cardboard



ORDERING SPECIFICATIONS

Item no.	Description	Pack	Quantity	Box
000 252 081	Neonate face mask (Size 1)	1		20
000 252 082	Infant face mask (Size 2)	1		20
000 252 083	Toddler face mask (Size 3)	1		20
000 252 084	Small adult/child face mask (Size 4)	1		20
000 252 085	Medium adult face mask (Size 5)	1		20
000 252 086	Large adult face mask (Size 6)	1		20



Ambu

Ambu Australia Pty Ltd
Unit 5, 2 Daydream Street
Warriewood, NSW 2102
Tel.: 1300 233 118
Fax: 1300 156 526
aussales@ambu.com
www.ambuaustralia.com.au