

Ambu® Specimen Suction Set



Key Benefits

- Female connector on inlet tube connects to catheter/bronchoscope and male connector fits to suction source.
- Transparent container provides easy visual examination of the specimen.
- Extended inlet length helps retain the specimen in the container.
- Container available in 3 sizes.
- Additional closure cap provided for safe transportation.

Ambu® Specimen Suction Set

For obtaining a tracheobronchial mucus specimen for microbiological examination.

For use with a separate suction catheter or flexible optic scope.

Ordering Specifications

Product Code	Description	Quantity per pouch	Quantity per bag	Box quantity
G05256	Specimen Suction Set 10 mls	1	50	150
G05502	Specimen Suction Set 25 mls	1	50	150
G05082	Specimen Suction Set 50 mls	1	50	150

Ambu Australia Pty Ltd
Customer Services
Tel.: 1300 233 118
Fax. 1300 156 526
aussales@ambu.com

Ambu Asia (Sales office)
Customer Services
Tel.: +604 252 9100
Fax. +604 644 8932
csasia@ambu.com

Ambu
Ideas that work for life