

Ambu® Mark III Adult Military



The original Ambu Bag with its unique double wall concept is designed for toxic atmospheres. The Ambu Mark III is a latex-free reusable resuscitator designed for manual ventilation of adults and children with a minimum body weight of approx. 15kg (3 years). The Mark III has the added advantage of being constructed with an outer made from butyl and an

optional gas filter adapter which is ideal for manual resuscitation in toxic atmospheres, e.g. fire and chemical accidents. The Ambu Mark III can be repeatedly autoclaved at 134°C (except for the face mask that is autoclaved at 121°C). It comes with a built-in pressure limitation of approx. 70cm H₂O for patient safety.

Main Technical Data

Mark III Adult Military

Max. tidal volume approx: 1300ml

Size (Length x diameter): 300 x 125mm (11.8" x 4.9")

Weight of resuscitator with patient valve but without reservoir and mask approx: 335g

Maximum frequency of ventilation: Depends on the inspiratory volume used. The rate of expansion of the bag is more than adequate for all frequencies recommended for practical use

Patient connector: 22/15mm (ISO)

Expiratory connector: 30mm male (ISO)

Main features and benefits

■ Excellent feeling of lung compliance

The unique double wall design provides a good sensation of ventilation pressure and lung compliance when the bag is compressed.

■ Free positioning of the bag

The patient swivel-valve enables free positioning of the bag without removing the mask from the patient.

■ Reduced fatigue during manual ventilation

The hand strap ensures a good grip, which helps to reduce fatigue during manual ventilation; also of benefit for users with small hands.

■ Easy to clean and sterilize

Military Mark III resuscitator is easy to disassemble. An effective procedure for cleaning and sterilization can be used. Can be autoclaved at 134°C.



ORDERING INFORMATION

Ambu Military Mark III resuscitator, complete with butyl outer cover, patient valve, gas filter adapter (DIN norm 3182), adult transparent face mask size 5 with black cuff and special hanging strap with snap hook

A233 004 000

Special military mask

ORDERING INFORMATION

Ambu Transparent Face Mask size 2 Child

Ambu Transparent Face Mask size 3/4 Child/Adult

Ambu Transparent Face Mask size 5 Adult

CAT. NO.

A000 016 000

A000 012 000

A000 017 000

